

# Starla



**Stannah**

**Arthur's always leaving his reading glasses around the house. But now that he's got a Stannah stairlift, it's easy for him to search high and low for them. Mind you, it would be even easier sometimes if he realised that they were on the top of his head all along.**



**We are intensely proud to have stayed a family-owned independent business. It has allowed us to remain true to our commitment to making the highest quality stairlifts.**

**Our Starla stairlift blends reliability with innovation. The intuitive and user-friendly design makes the Starla feel like it's built just for you.**

#### **The chair**

The intuitive, easy-to-use controls can be placed on the left or right arm to suit you and make operating the chair second nature. Sitting in a comfortable seat that supports your body, you simply move the control left or right to go up or down the stairs. When not in use the chair folds up effortlessly; simply raise the arms, seat and footrest to leave your stairs and hallway clear and tidy.

#### **The rail**

The Starla chair is available for straight or curved staircases. For straight stairs the Starla is fitted to a strong rail which is available in a silver or bronze finish. Our curved stairlifts are tailor-made to fit your individual staircase and ensure the best fit possible. The twin rails are available in a pearl colour as standard, with an additional five colour options.

Our rails are fitted to the stairs, not to the wall and installation is quick and clean, meaning minimum fuss and disruption in your home.

#### **The service**

The reliability of our stairlifts comes from the quality of the engineering and design. For your added peace of mind, all our new and reconditioned products are backed up by a full warranty and service package. Just like any piece of machinery, we recommend that your stairlift is serviced annually by a trained engineer. We will arrange with you to come and visit to give the stairlift a thorough and complete service, putting it through its paces.

You won't have to wait for service or help as our nationwide teams of engineers are locally-based. They are on standby 24-hours a day, 365 days a year to respond quickly in the unlikely event that you might need assistance.

### Upholstery material

Choose an all-fabric or all-vinyl chair or add a dark or light wood trim.

### Self-charging battery

The Starla always has the optimum amount of power, so you get a smooth ride with a gentle start and stop. It even works if there is a power cut.

### Vinyl upholstery colours



### Woven upholstery colours



### Swivel seat

Using the levers enables you to swivel the seat easily, making it simple to get off the chair. Alternatively, you can choose to have a powered swivel seat, which turns the chair to face the landing automatically.

### Optional two-way powered swivel

If space is limited at the bottom of the stairs, the two-way powered swivel will enable you to get off the chair with ease.

### Powered footrest button

This button helps with easy folding and is ideal if you have difficulty lifting your feet when getting onto the stairlift.

### Easy-to-use seat belt

Fastening the seat belt with one hand is simple and easy and gives you added security.

### Straight rails

Fitted to the stairs and not to the wall, the straight rail is available in a silver or bronze finish.

### Curved rails

Fitted to the stairs and not to the wall, the twin curved rails are available in a pearl colour as standard with an additional five colour options.



### Easy-to-use controls

Our new, easier-to-use, slimline control makes operating your Starla effortless.

### Arm isolation

For your safety, your Starla will not operate using the arm control until the armrest nearest the bottom of the staircase is fully lowered.

### Seat-load sensor

The extra safe sensor feature means that your stairlift won't move until you're sitting safely in the seat.

### Wall control

You can use the wall controls to call the stairlift for yourself or send it to someone else.

### Safety edges

If there is something blocking the chair on the stairs, the safety edges recognise this and the lift will stop automatically.



Curved model featured





1. Your Starla folds flat against your wall when not in use, giving you and others plenty of space.



2. Simply unfold the armrest and the seat. Lower the footrest by pressing the powered footrest button and your Starla is ready to use. Fasten the seatbelt for added security.



3. Gently move the slimline control and your Starla will glide smoothly up or down your stairs.



4. At the top of your stairs, you can swivel your Starla chair around to get off safely.



5. When not in use, you can leave your Starla in the swivelled position to act as a safety barrier.



6. Use the wall-mounted remote control to call your Starla to the top or bottom of your stairs.

